

ESWT Transmitters for R-SW Handpiece »ultra«

STORZ MEDICAL



PERI-ACTOR® Transmitters (according to Swart/Di Maio) Treatment of Fascia (max. pressure up to 3 bar)

PERI-ACTOR® »knuckle«



Treatment of all medium and big fasciae e.g.

- Back fascia, Plantar fascia
- Trigger bands
- Fasciae in the shoulder and neck area

PERI-ACTOR® »scoop«



Treatment on or below the bone edges e.g.

- Pelvis from the front
- Tibial crest
- Elbow

PERI-ACTOR® »sphere«



Punctual treatment of the fasciae e.g.

- Plantar fascia
- Arch of the foot
- Back of the hand
- Trigger points
- Cervical muscles

PERI-ACTOR® »scraper«



Treatment of all fasciae of the locomotor system

- Crural fascia
- Fascia thoracolumbalis
- Trigger bands
- Fascia lata

SPINE-ACTOR® Transmitters

Treatment of Paravertebral Muscles (max. pressure up to 3 bar)

SPINE-ACTOR® 1



- Span width: 20 mm

SPINE-ACTOR® 2



- Span width: 27.5 mm

SPINE-ACTOR® 3



- Span width: 35 mm

Treatment of the vertebral joints of the entire spinal column in the cervical, thoracic and lumbar areas e.g. treatment of the paravertebral muscles, dysfunctions in the entire spinal column, muscular imbalances and pain conditions as well as trigger points